

AChE check Control high/low

Professional quality management made easy

www.securetec.net



AChE check Control

AChE check Control makes the results of AChE testing more reliable and conclusive. The use of controls is essential for compliance with legal requirements and standards.

- ✓ Ease of use
- ✓ High similarity to patient samples
- ✓ Predefined targets and tolerance ranges for use with ChE check mobile
- ✓ Can also be used with AChE tests from other manufacturers

AChE check Control consists of two lyophilized whole blood control samples. One control has an acetylcholinesterase activity within the normal range. The other control sample is adjusted to a lower/inhibited activity.

This ensures that the medically relevant decision ranges are covered.



Item number: CH1300



Securetec Detektions-Systeme AG

Lilienthalstraße 7 · 85579 Neubiberg · Germany

T +49 89 203080-1651 · F +49 89 203080-1652

info@securetec.net · www.securetec.net

© Securetec Detektions-Systeme AG 2016 · 70531-EN-v01-2016-06-25 · Image source: fotolia